Thyroglobulin

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Thyroglobulin

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of thyroglobulin in serum or heparinized plasma, as an aid in monitoring patients who have undergone thyroidectomy.

Catalog Number: **L2KTY2** (200 tests)

Test Code: **TG** Color: **Dark Gray**

The presence of autoantibodies against thyroglobulin (anti-TG) can interfere with the determination of thyroglobulin in serum. Anti-TG antibodies can be quantitated by using the IMMULITE 2000 Anti-TG Ab kit (L2KTG).

When monitoring patients over time, TG values obtained with different methods cannot be used interchangeably due to the differences in methods and reagent specificity. The laboratory should provide the physician with the method identity along with the result. Serial serum TG measurements should be made with the same method.

Summary and Explanation

Thyroglobulin (TG) is a heterogeneous iodoglycoprotein which has a molecular mass of approximately 660,000 daltons. Thyroglobulin is normally synthesized in the follicular cells of the thyroid gland, under the influence of thyrotropin, and represents the precursor to thyroxine and the other iodothyronines.^{2,12}

The expected upper limit of normal for circulating thyroglobulin is approximately 40 to 60 ng/mL, with a median of 5 to 10 ng/mL.^{5,11} Somewhat higher values are encountered in newborns^{7,10} and during the third trimester of pregnancy. Thyroglobulin levels also tend to be elevated in regions of endemic goiter.

The major clinical applications for measurement of this prohormone derive from the fact that functioning thyroid tissue appears to be the only source of circulating thyroglobulin. Accordingly, thyroglobulin determinations have been widely used to complement radioiodine scanning and other techniques (such as ultrasound or immunohistochemical staining) as an aid in identifying the presence or absence of functioning thyroid tissue, or an increase in such tissue relative to an individually established baseline. The differential diagnosis of congenital hypothyroidism is a well-established context of use for this application of serum thyroglobulin measurement.

Congenital Hypothyroidism:

Thyroglobulin determinations have been used, sometimes in conjunction with ultrasound and radioiodine scanning, to help clarify the type of thyroid defect in previously diagnosed congenital hypothyroidism.^{3,4,6,7,10,13} Very low or undetectable thyroglobulin levels are expected in infants born without thyroid tissue (thyroid agenesis), whereas higher, but widely varying levels are generally encountered in infants with hypoplastic thyroid glands, ectopic thyroid tissue, dysmorphogenic goiter, congenital TBG deficiency or transient hypothyroidism.

Other Applications: Thyroglobulin measurements may also be of value in helping to distinguish subacute thyroiditis from thyrotoxicosis caused by covert administration of thyroid hormones. In the latter event, low levels of thyroglobulin are expected due to thyroid hormone suppression of thyrotropin.¹¹

Principle of the Procedure

IMMULITE 2000 Thyroglobulin is a solid-phase, chemiluminescent immunometric assay.

Incubation Cycles: 1 × 60 minutes

Specimen Collection

EDTA: Should not be used as an anticoagulant.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Thyroglobulin has not been tested with all possible variations of tube types.

Volume Required: 50 µL serum or heparinized plasma

Storage: 3 days at 2–8°C or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹⁸⁻²⁰

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H302 + H312, H412

P280, P273, P301 + P312, P302 + P312, P501

Warning! Harmful if swallowed or in contact with skin. Harmful to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations. **Contains:** sodium azide; Thyroglobulin Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. The barcode labels are needed for the assay.

Thyroglobulin Bead Pack (L2TY12)

With barcode. 200 beads, coated with anti-ligand. Stable at 2–8°C until expiration date.

L2KTY2: 1 pack

Thyroglobulin Reagent Wedge (L2TYA2)

With barcode. 11.5 mL of ligand-labeled anti-thyroglobulin murine monoclonal antibody and alkaline phosphatase (bovine calf intestine) conjugated to sheep polyclonal anti-thyroglobulin antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KTY2: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Thyroglobulin Adjustors (LTYL, LTYH)

Two vials (Low and High) of lyophilized thyroglobulin in a nonhuman serum matrix, with preservative. At least 30 minutes before use: Reconstitute each vial with **4.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Aliquot and freeze. Stable at –20°C for 2 months after reconstitution.

L2KTY2: 1 set

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Thyroglobulin Sample Diluent (L2TYZ):

For the on-board dilution of high samples. 25 mL of concentrated (ready-to-use) thyroglobulin-free nonhuman serum matrix. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2TYZ: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LTYCM: Tri-level, nonhuman serum-based control module

Also required

Sample transfer pipets, distilled or deionized water, controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:

2 weeks

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of thyroglobulin.

Expected Values

One hundred and ten serum samples from apparently healthy subjects were tested using the IMMULITE 2000 Thyroglobulin procedure. The results ranged from 0.73 to 84 ng/mL, with a median of 11 ng/mL and a 2.5th and 97.5th percentiles of 1.6 and 59.9 ng/mL, respectively.

Thirty serum samples from patients with Hashimoto's disease were tested using IMMULITE 2000 Thyroglobulin. The thyroglobulin values ranged from nondetectable to >300 ng/mL. Five of the 30 samples (17%) showed thyroglobulin levels greater than 55 ng/mL.

Another 30 serum samples from patients with Graves' disease were tested using IMMULITE 2000 Thyroglobulin. The thyroglobulin values ranged from nondetectable to >300 ng/mL. Three of the 30 samples (10%) manifested thyroglobulin levels greater than 55 ng/mL.

In patients who have undergone thyroidectomy, thyroglobulin concentrations should theoretically approach zero or approximate the sensitivity of the thyroglobulins assay. Any increase in the thyroglobulin level of such patients should be further interpreted in conjunction with other available diagnostic test results and the patient's clinical presentation.

Based on the studies of healthy subjects, the IMMULITE 2000 Thyroglobulin's cutoff was determined to be 55 ng/mL. Specimens with thyroglobulin levels greater than 55 ng/mL are considered elevated and require further testing regarding the thyroid status of the patient.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Because thyroglobulin autoantibodies (TgAb) can interfere with both competitive immunoassays and immunometric assays for thyroglobulin, all patients should be screened for TgAb by a sensitive immunoassay: recovery studies are not adequate for ruling out interference by TgAb.¹⁴⁻¹⁶ If the laboratory chooses to report thyroglobulin results in TgAb-positive sera, each report should include a warning to the physician to interpret the thyroglobulin result with caution, as the presence of TgAb may cause underestimation of the thyroglobulin level as measured by immunometric assays.¹⁵ Thyroglobulin antibodies (TgAb) are present in the majority of patients with Hashimoto's thyroiditis, but also in approximately 5% of healthy individuals.⁵

At least 6 weeks should elapse after thyroidectomy or ¹³¹I treatment.¹¹ Some reports have indicated that thyroglobulin levels may remain elevated for several months following successful treatment. In this case, serial determinations assessed relative to a post-treatment baseline established for the patient may still be of value in monitoring.

It is important to obtain thyroglobulin measurements before administering ¹³¹I for scanning, since the scan may result in the release of large amounts of thyroglobulin from even a very small number of cells. Likewise, samples for thyroglobulin assay should not be collected prematurely following procedures such as needle biopsy that are likely to cause a transient elevation in the thyroglobulin level.¹¹

Although thyroglobulin measurements are of established value in the differential diagnosis of congenital hypothyroidism, they are *not* recommended for use in screening for this condition since a wide spectrum of thyroglobulin levels are obtained for this condition.

EDTA Plasma has an effect on the measurement of thyroglobulin in the procedure.

Immunologically inactive thyroglobulin may not be determined by this assay.

In patients receiving thyroid suppression therapy, thyroglobulin may not be reliably detected.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL.

Conversion Factor:

ng/mL $\times$ 1.515 $\rightarrow$ pmol/L
(CRM 457; Molecular mass 660 kDa)

Calibration Range: Up to 300 ng/mL
(455 pmol/L)

Standardized against Certified Reference Material for human thyroglobulin (CRM 457) in terms of the Community Bureau of Reference of the European Commission.

Analytical Sensitivity: 0.2 ng/mL
(0.3 pmol/L)

Functional Sensitivity: 0.9 ng/mL
(1.4 pmol/L)

High-dose Hook Effect:

None up to 500,000 ng/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three thyroglobulin solutions (380, 840 and 1,680 ng/mL) were assayed. (See "Recovery" table for representative data.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 μ L/mL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Specificity: The assay is specific for thyroglobulin, with no detectable crossreactivity to alpha-fetoprotein (AFP), FSH, TSH, diiodothyronine (T2), triiodothyronine (T3) and thyroxine (T4). (See "Specificity" table.)

Method Comparison: The assay was compared to a commercially available indirect solid phase enzyme immunometric assay for thyroglobulin (Kit A) on 150 samples from patients with Hashimoto's disease, Graves disease, patients who had undergone thyroidectomy, as well as euthyroid individuals. The thyroglobulin concentrations of these patients covered the entire calibration range of the assay. The data are tabulated in reference to the respective assays' suggested cutoffs.

IMMULITE 2000 Thyroglobulin				
Kit A	>55 ng/mL	$\leq$ 55 ng/mL	Relative Sensitivity	Relative Specificity
>50 ng/mL	10	0	100%	92.9%
$\leq$ 50 ng/mL	10	130		

Agreement: 93.3%

95% Confidence Limits for Relative Sensitivity and Specificity, respectively: 69.2% – 100% and 87.3% – 96.5%,

Clinical Performance: A clinical study at a university medical center in the midwestern United States included 56 thyroid disease patients (41 females and 15 males, with an age range from 20 to 72 years) who provided from 2 to 6 serial specimens in the course of their treatments that included partial, subtotal or total thyroidectomy, and other medical interventions. The results of this study showed that IMMULITE 2000 Thyroglobulin measurements were consistent with the thyroid disease conditions of these patients.

References

- 1) Bodlaender P, et al. Sensitive radioimmunological screening test for anti-thyroglobulin autoantibodies. Clin Chem 1978;24:272-4.
- 2) Cody V. Thyroglobulin and thyroid hormone synthesis. Endocr Res 1984;10:73-88.
- 3) Czernichow P, et al. Plasma thyroglobulin measurements help determine the type of thyroid defect in congenital hypothyroidism. J Clin Endocrinol Metab 1983;56:242-5.
- 4) Dammacco F, et al. Serum thyroglobulin and thyroid ultrasound studies in infants with congenital hypothyroidism. J Pediatr 1985;106:451-3.
- 5) Feldt-Rasmussen U. Serum thyroglobulin and thyroglobulin autoantibodies in thyroid diseases. Allergy 1983;38:369-87.
- 6) Gons MH, et al. Concentration of plasma thyroglobulin and urinary excretion of iodinated material in the diagnosis of thyroid disorders in congenital hypothyroidism. Acta Endocrinol 1983;104:27-34.
- 7) Ket JL, et al. Serum thyroglobulin levels: physiological decrease in infancy and the absence in athyroidism. J Clin Endocrinol Metab 1981;58:1301-3.
- 8) Lo Gerfo P, et al. Thyroglobulin in benign and malignant thyroid disease. JAMA 1979;241:923-5.
- 9) Pacini F, et al. Serum thyroglobulin in thyroid carcinoma and other thyroid disorders. J Endocrinol Invest 1980;3:283-92.
- 10) Pinchera A, et al. Recent studies on the application of serum thyroglobulin measurement in thyroid disease. In: Eggo M, Burrow GN, editors. Thyroglobulin – the prothyroid hormone. New York: Raven Press, 1985: 307-16.
- 11) Refetoff S, Lever EG. The value of serum thyroglobulin measurement in clinical practice. JAMA 1983;250:2352-57.
- 12) Van Herle AJ, Vassart G, Dumont JE. Control of thyroglobulin synthesis and secretion. N Engl J Med 1979;301:239-49 and 307-14.
- 13) Walfish PG. Thyroid function in pediatrics. In: Hicks J, Boeckx RL, editors. Pediatric clinical chemistry, Philadelphia: W.B. Saunders, 1984: 170-239.
- 14) Spencer CA, Wang C-C. Thyroglobulin measurement: techniques, clinical benefits, and pitfalls. Endocrinol Metab Clin North Am 1995;24: 841-63.
- 15) Spencer CA, Takeuchi M, Kazarosyan M. Current status and performance

goals for serum thyroglobulin assays. Clin Chem 1996;42:164-73. 16) Spencer CA. Recoveries cannot be used to authenticate thyroglobulin (Tg) measurements when sera contain Tg autoantibodies. Clin Chem 1996;42:661-3. 17) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 18) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8. 19) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 20) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	10	0.48	4.8%	0.56	5.6%
2	59	3.5	5.9%	3.6	6.1%
3	13	0.82	6.3%	0.85	6.5%
4	20	1.0	5.0%	2.0	10%
5	279	19	6.8%	20	7.2%

Linearity (ng/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	16 in 16	63	—	—
	8 in 16	31	32	97%
	4 in 16	15	16	94%
	2 in 16	7.3	7.9	92%
	1 in 16	3.7	3.9	95%
2	16 in 16	82	—	—
	8 in 16	37	41	90%
	4 in 16	19	21	91%
	2 in 16	10	10	100%
	1 in 16	4.7	5.1	92%
3	16 in 16	96	—	—
	8 in 16	44	48	92%
	4 in 16	21	24	88%
	2 in 16	10	12	83%
	1 in 16	5.1	6.0	85%
4	16 in 16	140	—	—
	8 in 16	70	70	100%
	4 in 16	34	35	97%
	2 in 16	17	18	94%
	1 in 16	8.3	8.8	94%

Specificity

Compound ¹	Added ² ng/mL	Apparent ng/mL ³	% Cross- reactivity ⁴
AFP	19,000	ND	ND
FSH	200	0.11	ND
TSH	560	0.09	ND
T2	48,000	ND	ND
T3	10,000	ND	ND
T4	10,000	ND	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.

© 2008 Siemens Healthcare Diagnostics. All rights reserved.

Made in: UK

Siemens Healthcare
Diagnostics Products Ltd.
Glyn Rhonwy, Llanberis,
Gwynedd LL55 4EL
United Kingdom

2015-04-22

PIL2KTY – 22

Changes in this Edition:

cc#EU22516: 1) Under Warnings and Precautions, added biohazard symbol and potential biohazard statement, caution statement, required GHS pictograms/signal words/hazard and precautionary codes and statements, offending chemicals and sources. 2) Under References, added citations from the potential biohazard statement. 3) Changed "Origin: UK" to "Made in: UK". 4) Updated the manufacturer address to include Glyn Rhonwy. 5) In Understanding the Symbols, revised the symbol for catalog number, added a GHS pictogram for health hazard, and updated existing symbols with GHS pictograms and terminology for exclamation mark, corrosion, skull and crossbones, and environment.

Understanding the Symbols

Understanding the Symbols	En English
Erklärung der Symbole	De Deutsch
Descripción de los símbolos	Es Español
Explication des symboles	Fr Français
Definizione dei simboli	It Italiano
Descrição dos símbolos	Pt Português

The following symbols may appear on the product labeling: / Die folgenden Symbole können auf dem Produktetikett erscheinen: / Los siguientes símbolos pueden aparecer en la etiqueta del producto: / Les symboles suivants peuvent apparaître sur les étiquettes des produits: / Sull'etichetta del prodotto possono essere presenti i seguenti simboli: / Os seguintes símbolos podem aparecer no rótulo dos produtos:

Symbol Definition

En: *In vitro* diagnostic medical device

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CONTROL AB En: Control Antibody