

# m-PIMA™ HIV-1/2 VL

Viral load monitoring allows patients to stay on more affordable, less complex first-line regimens, with clear benefits for both patient and the healthcare system<sup>1</sup>.

Challenges for continuous HIV Viral Load monitoring include access, timeliness of result, and retention

- **Access:** Cost and complexity of currently available technologies limit availability in resource-poor settings<sup>2</sup>.
- **Long turnaround time:** Centralized testing causes lengthy delays in receiving results and disruptions in the treatment cascade<sup>2,3</sup>.
- **Loss to follow-up:** <50% of virally suppressed patients receive a second, or follow-up HIV Viral Load test<sup>4</sup>.



# m-PIMA™ HIV-1/2 VL

## ACTIONABLE RESULTS IN LESS THAN 70 MINUTES WHILE THE PATIENT IS STILL PRESENT

- **Extending VL Reach:** Early identification of treatment failure/lack of adherence can prevent transmissions, halt disease progression and reduce mortality.
- **Actionable Results at the Point of Care:** Results in <70 minutes while the patient is still present.
- **Potential Reduction of Patients Lost to Follow Up:** Improving the standard of care by bringing care to where it is needed the most (rural health centers).
- **Manage decentralized programs:** m-PIMA™ has connectivity enabled via CONNECT Universal Gateway.

**Product Specifications:** WHO recommends that all patients on ART receive a viral load test at 6 months and 12 months, and annually thereafter if the patient is stable on ART <sup>1</sup>, but very few patients receive that level of care.



CURRENT  
ALTERNATIVES  
**WEEKS  
OR MONTHS**

**VS**



POINT OF CARE  
**LESS THAN  
70 MINUTES**

The m-PIMA™ HIV-1/2 VL test is a quantitative nucleic acid amplification test designed for the quantification of HIV-1 groups M/N, and O and HIV-2 RNA in human plasma specimens.

### Accuracy

HIV-1 (group M/N and O) and HIV-2  
Sensitivity: ≥99%\*  
Specificity: ≥98%\*

### LOD (Limit of Detection)

HIV-1 group M: 342 copies/mL  
HIV-1 group O: 228 copies/mL  
HIV-2: 364 copies/mL

### LLOQ (Lower Limit of Quantitation)

800 copies/mL

\*at 1,000 copies/mL threshold for virologic failure

### Time to Results

<70 minutes

### Easy to perform

Sample volume: 50 µL EDTA plasma generated from venous or capillary whole blood.



## VIRAL LOAD MONITORING HAS THE POTENTIAL TO SUPPORT THE REALIZATION OF THE 95-95-95 GOALS.

### ORDER INFO

PART	ORDER NUMBER	QTY/BOX
m-PIMA™ Analyser	27030R001	1
m-PIMA™ HIV-1/2 VL	270150050	50
USB Printer	27040R007	1
Connect Universal Gateway	ACEU01	1

NOT CLEARED/APPROVED OR FOR SALE IN USA.

1. Bonner et al. Viral Load Monitoring as a Tool to Reinforce Adherence: A Systematic Review. Journal of Acquired Immune Deficiency Syndromes. (Sep 2013) 64(1):74-78.
2. Point-of-care viral load tests to detect high HIV viral load in people living with HIV/AIDS attending health facilities. Cochrane Database of Systematic Reviews (2022).
3. Pham, M. D., Nguyen, H. V., Anderson, D., Crowe, S., & Luchters, S. Viral load monitoring for people living with HIV in the era of test and treat: progress made and challenges ahead – a systematic review. BMC Public Health. (June 2022); 22(1):1203.
4. WHO. Toolkit: HIV molecular diagnostics toolkit to improve access to viral load testing and infant diagnosis: HIV treatment and care (2019). HIV molecular diagnostics toolkit to improve access to viral load testing and infant diagnosis. Geneva: World Health Organization; 2019. Licence: CC BY-NC-SA 3.0 IGO.

© 2025 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. RDx-24026950-3.0 03/25

